

Anti-Bcl2L14/Bclg Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Bzr8

Item Number: CM0003616

Introduction

Rabbit monoclonal antibody against human BCL2L14/BCLG, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for apoptosis research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	BCL2L14/BCLG Rabbit mAb
Remarks/Alias	BCLG; BCL2L14/BCLG
Species	Human
GeneID (Human)	79370
GeneID	79370
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human BCL2L14/BCLG (Q9BZR8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZR8
Protein Weight	37kDa
Shipping	Ice bag