

Anti-Bcl2L14/Bclg Rabbit Polyclonal Antibody For Wb, If, Elisa - Q9Bzr8

Item Number: CM0028050

Introduction

Rabbit polyclonal antibody against human BCL2L14/BCLG (apoptosis facilitator Bcl-2-like 14). Validated for WB, IF/ICC, ELISA in human, mouse, rat. Suitable for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	BCL2L14/BCLG Rabbit pAb
Remarks/Alias	BCLG; BCL2L14/BCLG
Species	Human
Gene ID (Human)	79370
Gene ID	79370
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-327 of human BCL2L14/BCLG (NP_620049.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZR8
Protein Weight	37kDa
Shipping	Ice bag