


Anti-Bcl2L15 Rabbit Monoclonal Antibody For Wb, Elisa - Q5Tbc7

Item Number: CM0006703



Introduction

Recombinant rabbit monoclonal antibody against human BCL2L15 (Bfk, C1orf178). Validated for WB and ELISA. Useful for apoptosis, alternative splicing, and proteomics research applications.

[Learn More](#)

Parameter	Value
Product Name	BCL2L15 Rabbit mAb
Remarks/Alias	Bfk; C1orf178; BCL2L15
Species	Human
GeneID (Human)	440603
GeneID	440603
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human BCL2L15 (NP_001010922.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5TBC7
Protein Weight	18kDa
Shipping	Ice bag