

Anti-Bcl3 Polyclonal Antibody For Wb, If/Icc, Elisa - P20749

Item Number: CM0014684



Introduction

Rabbit polyclonal antibody against human BCL3, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat BCL3. Suitable for studies of NF-kappa-B regulation and oncogenic signaling.

[Learn More](#)

Parameter	Value
Product Name	BCL3 Rabbit pAb
Remarks/Alias	BCL4; D19S37; BCL3
Species	Human
GeneID (Human)	602
GeneID	602
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-400 of human BCL3 (NP_005169.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20749
Protein Weight	48kDa
Shipping	Ice bag