

Anti-Bcl6 Monoclonal Antibody () For Intracellular Flow Cytometry - P41182

Item Number: CM0012275



Introduction

Rabbit anti-human BCL6 monoclonal antibody directly conjugated to for intracellular flow cytometry. Validated for FC (intra), targets a transcriptional repressor critical in germinal center biology and lymphoma research. Human cross-reactivity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human BCL6 mAb
Remarks/Alias	BCL6; BCL5; BCL6A; LAZ3; ZBTB27; ZNF51; B-cell CLL/lymphoma 6
Species	Human
GeneID	604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 235-350 of human BCL6 (NP_001124317.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P41182
Protein Weight	72kDa/78kDa
Shipping	Ice bag