

Anti-Bcl6 Rabbit Monoclonal Antibody For Flow Cytometry And Elisa - P41182

Item Number: CM0006651



Introduction

Recombinant anti-human BCL6 rabbit monoclonal antibody validated for intracellular flow cytometry and ELISA. Ideal for studying B-cell lymphoma 6 protein function and localization.

[Learn More](#)

Parameter	Value
Product Name	BCL6 Rabbit mAb
Remarks/Alias	BCL5; LAZ3; BCL6A; ZNF51; ZBTB27; BCL6
Species	Human
GeneID (Human)	604
GeneID	604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 235-350 of human BCL6 (NP_001124317.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra), ELISA
Cross Reactivity	Human
SWISS	P41182
Protein Weight	79kDa
Shipping	Ice bag