


Anti-Bcl9 Rabbit Polyclonal Antibody For Wb And Elisa - 000512

Item Number: CM0022392



Introduction

Rabbit polyclonal antibody against human Bcl9, a Wnt pathway proto-oncogene. Validated for Western blot (WB) and ELISA, with cross-reactivity to mouse. Immunogen: N-terminal recombinant protein. Ideal for cancer and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Bcl9 Rabbit pAb
Remarks/Alias	LGS; Bcl9
Species	Human
GeneID (Human)	607
GeneID	607
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human Bcl9. (NP_004317.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O00512
Protein Weight	149kDa
Shipping	Ice bag