

Anti-Bcl9 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9D219

Item Number: CM0028989

Introduction

Anti-Bcl9 rabbit polyclonal antibody suitable for WB, IHC-P, ELISA. Recognizes human, mouse, rat Bcl9, a ~149kDa Wnt pathway nuclear coactivator. Supports cancer and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Bcl9 Rabbit pAb
Remarks/Alias	Gm130; 2610202E01Rik; 8030475K17Rik; A330041G23Rik; Bcl9
Species	Mouse
GeneID (Human)	607
GeneID	77578
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1278-1377 of mouse Bcl9 (NP_084209.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9D219
Protein Weight	149kDa
Shipping	Ice bag