

Anti-Bcl9 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9D219

Item Number: CM0023964

Introduction

Rabbit polyclonal antibody against mouse B-cell CLL/lymphoma 9 (Bcl9) protein. Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for Wnt signaling research.

[Learn More](#)

Parameter	Description
Product Name	Bcl9 Rabbit pAb
Remarks/Alias	Gm130; 2610202E01Rik; 8030475K17Rik; A330041G23Rik; Bcl9
Species	Mouse
Gene ID (Human)	607
Gene ID	77578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-500 of mouse Bcl9 (Q9D219)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9D219
Protein Weight	149kDa
Shipping	Ice bag