


Anti-Bclaf1 Polyclonal Antibody For Wb And Elisa - Q9Nyf8

Item Number: CM0024202



Introduction

High-quality rabbit polyclonal antibody against BCLAF1 (BTF), suitable for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat BCLAF1. Ideal for apoptosis and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	BCLAF1 Rabbit pAb
Remarks/Alias	BTF, bK211L9.1, BCLAF1
Species	Human
Gene ID (Human)	9774
Gene ID	9774
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 330-552 of human BCLAF1 (NP_001070909.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYF8
Protein Weight	106kDa
Shipping	Ice bag