

Anti-Bcmo1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hay6

Item Number: CM0018310

Introduction

BCMO1 Rabbit Polyclonal Antibody for WB and ELISA applications. Validated in human, mouse, and rat samples. Targets beta-carotene dioxygenase 1, a key enzyme in retinal production and vision.

[Learn More](#)

Parameter	Value
Product Name	BCMO1 Rabbit pAb
Remarks/Alias	BCO; BCDO; BCMO; BCDO1; BCMO1
Species	Human
Gene ID (Human)	53630
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 268-547 of human BCMO1 (NP_059125.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HAY6
Protein Weight	63kDa
Shipping	Ice bag