


Anti-Bcr Polyclonal Antibody For Wb, Ihc-P, Elisa - P11274

Item Number: CM0015114



Introduction

High-quality anti-BCR rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat BCR. Ideal for studying small GTPase signaling and synaptic biology.

[Learn More](#)

Parameter	Value
Product Name	BCR Rabbit pAb
Remarks/Alias	ALL; CML; PHL; BCR1; D22S11; D22S662; BCR
Species	Human
GeneID (Human)	613
GeneID	613
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11274
Protein Weight	143kDa
Shipping	Ice bag