

Anti-Bdkrb2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P30411

Item Number: CM0029619



Introduction

Rabbit polyclonal antibody against human BDKRB2 (B2 bradykinin receptor), validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat. Suitable for GPCR signaling research and diagnostic development.

[Learn More](#)

□	□
Product Name	BDKRB2 Rabbit pAb
Remarks/Alias	B2R; BK2; BK-2; BKR2; BRB2; BDKRB2
Species	Human
GeneID (human)	624
GeneID	624
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 302-391 of human BDKRB2 (NP_000614.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30411
Protein Weight	44kDa
Shipping	Ice bag