


Anti-Bdnf Monoclonal Antibody For If-P, Elisa - P23560

Item Number: CM0005654



Introduction

Rabbit monoclonal antibody against human BDNF, validated for IF-P and ELISA. Detects mouse and rat BDNF. Ideal for neuroscience research on neurotrophic signaling, synaptic plasticity, and neuronal survival.

[Learn More](#)

Parameter	Value
Product Name	BDNF Rabbit mAb
Remarks/Alias	ANON2; BULN2; BDNF
Species	Human
Gene ID (Human)	627
Gene ID	627
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 100-200 of human BDNF (P23560).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P23560
Protein Weight	28kDa
Shipping	Ice bag