


Anti-Bdnf Monoclonal Antibody (Mouse) For Wb, Elisa - P23560

Item Number: CM0003884



Introduction

High-quality mouse monoclonal antibody against human BDNF (P23560), validated for Western blot and ELISA. Detects mouse and rat BDNF. Ideal for neurobiology and synaptic plasticity studies.

[Learn More](#)

Parameter	Value
Product Name	BDNF Mouse mAb
Remarks/Alias	ANON2; BULN2; BDNF
Species	Human
GeneID (Human)	627
GeneID	627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-247 of human BDNF (P23560)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P23560
Protein Weight	28kDa
Shipping	Ice bag