


Anti-Bdnf Rabbit Polyclonal Antibody For Wb, Elisa - P23560

Item Number: CM0022466



Introduction

Rabbit polyclonal antibody targeting human BDNF (Brain-derived neurotrophic factor), suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience and neurotrophin research.

[Learn More](#)

Parameter	Value
Product Name	BDNF Rabbit pAb
Remarks/Alias	ANON2; BULN2; BDNF
Species	Human
GeneID (human)	627
GeneID	627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-247 of human BDNF (NP_733928.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23560
Protein Weight	28kDa
Shipping	Ice bag