

Anti-Beclin 1 Rabbit Monoclonal Antibody (Ko Validated) For Wb And Elisa - Q14457

Item Number: CM0005889



Introduction

KO validated rabbit monoclonal antibody against human Beclin 1, validated for WB and ELISA. Recognizes human, mouse, and rat Beclin 1. Suitable for autophagy and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Beclin 1 Rabbit mAb
Remarks/Alias	ATG6; VPS30; beclin1; Beclin 1; BECN1; beclin-1; 1
Species	Human
GeneID (human)	8678
GeneID	8678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human Beclin 1 (NP_003757.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14457
Protein Weight	52kDa
Shipping	Ice bag