


Anti-Best1 Polyclonal Antibody For Wb, Elisa - 076090

Item Number: CM0013493



Introduction

Rabbit polyclonal antibody targeting human BEST1 (Bestrophin-1), validated for WB and ELISA. Reacts with human and mouse samples. Ideal for retinal pigment epithelium and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	BEST1 Rabbit pAb
Remarks/Alias	ARB; BMD; BEST; RP50; VMD2; TU15B; Best1V1Delta2; BEST1
Species	Human
GeneID (Human)	7439
GeneID	7439
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O76090
Protein Weight	68kDa
Shipping	Ice bag