


Anti-Best1 Rabbit Polyclonal Antibody For Wb, Elisa - 076090

Item Number: CM0030348



Introduction

Rabbit polyclonal antibody against human Bestrophin-1 (BEST1). Validated for WB and ELISA; cross-reacts with human, mouse, rat. Ideal for retinal pigment epithelium research and anion channel studies.

[Learn More](#)

Parameter	Value
Product Name	BEST1 Rabbit pAb
Remarks/Alias	ARB; BMD; BEST; RP50; VMD2; TU15B; Best1V1Delta2; BEST1
Species	Human
GeneID (Human)	7439
GeneID	7439
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 171-270 of human BEST1 (NP_004174.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	076090
Protein Weight	68kDa
Shipping	Ice bag