

Anti-Bet1 Rabbit Polyclonal Antibody For Wb And Elisa - O15155

Item Number: CM0015472



Introduction

Rabbit polyclonal antibody against human BET1 (HBET1), validated in WB and ELISA. BET1 is a SNARE protein required for ER-to-Golgi vesicular transport, localized at the ER and cis-Golgi membranes. Ideal for studying intracellular trafficking.

[Learn More](#)

Parameter	Value
Product Name	BET1 Rabbit pAb
Remarks/Alias	HBET1; BET1
Species	Human
GeneID (Human)	10282
GeneID	10282
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-94 of human BET1 (NP_005859.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15155
Protein Weight	13kDa
Shipping	Ice bag