

Anti-Bet1L Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nym9

Item Number: CM0018642

Introduction

Rabbit polyclonal antibody targeting human BET1L (GS15), a Golgi SNARE protein. Suitable for WB and ELISA applications. BET1L is involved in retrograde vesicle transport and Golgi complex integrity. Ideal for research on Golgi trafficking and SNARE biology.

[Learn More](#)

Parameter	Value
Product Name	BET1L Rabbit pAb
Remarks/Alias	GS15; BET1L1; GOLIM3; HSPC197; BET1L
Species	Human
Gene ID (Human)	51272
Gene ID	51272
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human BET1L (NP_057610.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NYM9
Protein Weight	12kDa
Shipping	Ice bag