

Anti-Beta 2 Microglobulin (B2M) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Fc, Elisa - P61769

Item Number: CM0007359



Introduction

Validated anti-beta 2 Microglobulin Rabbit mAb for WB, IHC-P, IF/ICC, FC, ELISA. Recognizes human B2M (13kDa). Ideal for MHC class I antigen presentation and amyloidosis studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] beta 2 Microglobulin Rabbit mAb
Remarks/Alias	B2M; IMD43; beta-2-microglobulin
Species	Human
Gene ID	567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-119 of human β 2 Microglobulin (NP_004039.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P61769
Protein Weight	13kDa
Shipping	Ice bag