


Anti-Beta 2 Microglobulin Rabbit Mab For Wb, If-P, Ihc-P, Elisa - P61769

Item Number: CM0003212



Introduction

Rabbit monoclonal antibody against human beta 2 Microglobulin (B2M), suitable for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. A key tool for MHC class I antigen presentation studies.

[Learn More](#)

Parameter	Value
Product Name	beta 2 Microglobulin Rabbit mAb
Remarks/Alias	IMD43; beta 2 Microglobulin
Species	Human
GeneID (Human)	567
GeneID	567
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61769
Protein Weight	14kDa
Shipping	Ice bag