

Anti-Beta 2 Microglobulin Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P61769

Item Number: CM0015321



Introduction

Rabbit polyclonal antibody targeting human beta 2 Microglobulin, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for MHC I and immune research.

[Learn More](#)

Parameter	Value
Product Name	beta 2 Microglobulin Rabbit pAb
Remarks/Alias	IMD43; beta 2 Microglobulin
Species	Human
Gene ID (Human)	567
Gene ID	567
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-119 of human beta 2 Microglobulin (NP_004039.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P,IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61769
Protein Weight	14kDa
Shipping	Ice bag