

Anti-Beta-Defensin 26 (Defb26) Rabbit Polyclonal Antibody For If/Icc, Elisa - Q32Zg6

Item Number: CM0013142

Introduction

Defb26 Rabbit pAb polyclonal antibody validated for IF/ICC, ELISA. Cross-reacts with human. Raised in rabbit against recombinant rat Defb26 protein. For antimicrobial defensin research.

[Learn More](#)

Parameter	Value
Product Name	Defb26 Rabbit pAb
Remarks/Alias	DEFB126; C20orf8; DEFB-26; DEFB26; HBD26; bA530N10.1; hBD-26; Defb26
Species	Rat
Gene ID (Human)	81623
Gene ID	641629
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-128 of rat Defb26 (NP_001032595.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q32ZG6
Protein Weight	14kDa
Shipping	Ice bag