

Anti-Beta-Galactosidase (Glb1) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P16278

Item Number: CM0007386



Introduction

High-quality rabbit monoclonal antibody against human Beta-galactosidase (GLB1), validated for WB, IHC-P, and ELISA. Suitable for detecting GLB1 in lysosomal and elastogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Beta-galactosidase (GLB1) Rabbit mAb
Remarks/Alias	GLB1; EBP; ELNR1; MPS4B; beta-galactosidase
Species	Human
Gene ID (Human)	2720
Immunogen	Recombinant protein of human Beta-galactosidase (GLB1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P16278
Protein Weight	60kDa/72kDa/76kDa
Shipping	Ice bag