

Anti-Beta-Galactosidase (Glb1) Rabbit Polyclonal Antibody For Wb, Elisa - P16278

Item Number: CM0016241



Introduction

Rabbit polyclonal antibody against human GLB1 (beta-galactosidase/elastic receptor 1), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for research on lysosomal storage disorders and elastogenesis.

[Learn More](#)

Parameter	Specification
Product Name	GLB1 Rabbit pAb
Remarks/Alias	EBP; ELNR1; MPS4B; GLB1
Species	Human
Gene ID	2720 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human GLB1 (NP_000395.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16278
Protein Weight	76kDa
Shipping	Ice bag