


Anti-Beta Glucuronidase (Gusb) Rabbit Monoclonal Antibody For Wb And Elisa - P08236

Item Number: CM0002563



Introduction

This rabbit monoclonal antibody targets human beta-glucuronidase (GUSB), a lysosomal enzyme that degrades dermatan and keratan sulfates. Ideal for WB and ELISA applications using human samples. Immunogen corresponds to amino acids 1-100.

[Learn More](#)

Parameter	Value
Product Name	beta glucuronidase (GUSB) Rabbit mAb
Remarks/Alias	BG; MPS7; beta glucuronidase (GUSB)
Species	Human
GeneID (Human)	2990
GeneID	2990
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human beta-glucuronidase (GUSB) (NP_000172.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P08236
Protein Weight	75kDa
Shipping	Ice bag