

Anti-Bfar Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nzs9

Item Number: CM0025562

Introduction

Rabbit polyclonal antibody targeting human BFAR (RNF47), validated for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat samples. Suitable for apoptosis and ER stress research.

[Learn More](#)

Parameter	Value
Product Name	BFAR Rabbit pAb
Remarks/Alias	BAR; RNF47; BFAR
Species	Human
GeneID	51283
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human BFAR (NP_057645.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZS9
Protein Weight	53kDa
Shipping	Ice bag