

Anti-Bgn Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P21810

Item Number: CM0028482

Introduction

Rabbit polyclonal antibody against human biglycan (BGN), suitable for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Ideal for extracellular matrix and connective tissue research.

[Learn More](#)

Parameter	Value
Product Name	BGN Rabbit pAb
Remarks/Alias	PGI; MRLS; DSPG1; PG-S1; SEMDX; SLRR1A; BGN
Species	Human
Gene ID (Human)	633
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-368 of human BGN (NP_001702.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P21810
Protein Weight	42kDa
Shipping	Ice bag