

Anti-Bhlhe40 Polyclonal Antibody For Wb, Ihc-P, Elisa - O14503

Item Number: CM0027088

Introduction

BHLHE40 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA in human, mouse, and rat. Targets the circadian transcriptional repressor DEC1/SHARP2. Ideal for biological rhythm and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	BHLHE40 Rabbit pAb
Remarks/Alias	DEC1; HLHB2; BHLHB2; Clast5; SHARP2; STRA13; Stra14; SHARP-2; BHLHE40
Species	Human
Gene ID (Human)	8553
Gene ID	8553
Immunogen	Recombinant protein or recombinant fusion protein containing a sequence corresponding to amino acids 121-300 of human BHLHE40 (NP_003661.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14503
Protein Weight	46kDa
Shipping	Ice bag