


Anti-Bhmt Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q93088

Item Number: CM0013713



Introduction

Rabbit polyclonal antibody targeting human BHMT (betaine-homocysteine S-methyltransferase 1). Validated for WB, IF-P, and ELISA applications, with cross-reactivity to mouse and rat. Ideal for homocysteine metabolism studies. Protein weight 45kDa.

[Learn More](#)

Parameter	Value
Product Name	BHMT Rabbit pAb
Remarks/Alias	BHMT1; HEL-S-61p; BHMT
Species	Human
GeneID (human)	635
GeneID	635
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 177-406 of human BHMT (NP_001704.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93088
Protein Weight	45kDa
Shipping	Ice bag