


Anti-Bicd2 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Td16

Item Number: CM0018401



Introduction

Rabbit polyclonal antibody against BICD2 (SMALED2) for Western blot and ELISA. Reacts with human, mouse, and rat. BICD2 is an adaptor protein for dynein motor complex, involved in Golgi-ER transport and nuclear positioning.

[Learn More](#)

Parameter	Specification
Product Name	BICD2 Rabbit pAb
Remarks/Alias	SMALED2; SMALED2A; SMALED2B; bA526D8.1; BICD2
Species	Human
GeneID (human)	23299
GeneID	23299
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 440-620 of human BICD2 (NP_001003800.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TD16
Protein Weight	94kDa
Shipping	Ice bag