


Anti-Bid Rabbit Mab (Kd Validated) For Wb, Elisa - P55957

Item Number: CM0007631



Introduction

KD validated rabbit monoclonal antibody against human Bid (UniProt P55957, 22 kDa). Suitable for Western blotting and ELISA; ideal for apoptosis and mitochondrial pathway studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Bid Rabbit mAb
Remarks/Alias	FP497; [KD Validated] Bid
Species	Human
Gene ID	637
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human Bid (NP_001187.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P55957
Protein Weight	22kDa
Shipping	Ice bag