


Anti-Bif-1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y371

Item Number: CM0016020



Introduction

Rabbit polyclonal antibody against human Bif-1 (Endophilin-B1, CGI-61), validated for WB, IHC-P, ELISA. Cross-reacts with mouse, rat. Bif-1 regulates mitochondrial dynamics, apoptosis, autophagy. For research.

[Learn More](#)

Parameters	Specification
Product Name	Bif-1 Rabbit pAb
Remarks/Alias	Bif-1; CGI-61; PPP1R70; dj612B15.2
Species	Human
GeneID (Human)	51100
GeneID	51100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-365 of human Bif-1 (NP_057093.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y371
Protein Weight	41kDa
Shipping	Ice bag