


Anti-Bif-1 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q9Y371

Item Number: CM0028409



Introduction

High-quality Bif-1 Rabbit Polyclonal Antibody validated for IF/ICC and ELISA. Targets human endophilin-B1 (SH3GLB1), a key regulator of mitochondrial dynamics, autophagy, and apoptosis. Ideal for IVD and research applications.

[Learn More](#)

Parameter	Value
Product Name	Bif-1 Rabbit pAb
Remarks/Alias	Bif-1; CGI-61; PPP1R70; dj612B15.2
Species	Human
Gene ID (human)	51100
Gene ID	51100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-365 of human Bif-1 (NP_057093.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y371
Protein Weight	41kDa
Shipping	Ice bag