


Anti-Bim Rabbit Monoclonal Antibody For Wb And Elisa - O43521

Item Number: CM0002482



Introduction

Validated Bim rabbit monoclonal antibody for WB and ELISA. Recognizes human, mouse, and rat Bim. Immunogen: synthetic peptide from N-terminal region. Ideal for apoptosis research. SWISS O43521.

[Learn More](#)

Parameter	Value
Product Name	Bim Rabbit mAb
Remarks/Alias	BAM; BIM; BOD; Bim
Species	Human
Gene ID (Human)	10018
Gene ID	10018
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Bim (O43521)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43521
Protein Weight	22kDa
Shipping	Ice bag