


Anti-Bim Rabbit Polyclonal Antibody For Ihc-P And Elisa - O43521

Item Number: CM0015044



Introduction

Rabbit polyclonal antibody against Bim (BCL2L11), a key apoptosis regulator, validated for IHC-P and ELISA in human samples. Ideal for cell death research.

[Learn More](#)

Parameter	Value
Product Name	Bim Rabbit pAb
Remarks/Alias	BAM; BIM; BOD; Bim
Species	Human
GeneID (Human)	10018
GeneID	10018
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-130 of human Bim (NP_619527.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	O43521
Protein Weight	22kDa
Shipping	Ice bag