


Anti-Bim Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43521

Item Number: CM0019470



Introduction

Rabbit polyclonal antibody against human Bim (BCL2L11) for apoptosis research. Validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying programmed cell death mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Bim Rabbit pAb
Remarks/Alias	BAM; BIM; BOD; Bim
Species	Human
GeneID (Human)	10018
GeneID	10018
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 100-198 of human Bim (NP_619527.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43521
Protein Weight	22kDa
Shipping	Ice bag