


Anti-Bin1 Monoclonal Antibody For Wb, Ihc-P, Elisa - 000499

Item Number: CM0003542



Introduction

BIN1 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA in Human, Mouse, Rat. Ideal for investigating membrane curvature, T-tubule formation, endocytosis, and MYC activity. Suitable for clinical diagnostic research and assay development.

[Learn More](#)

Parameter	Value
Product Name	BIN1 Rabbit mAb
Remarks/Alias	CNM2; AMPH2; AMPHL; SH3P9; BIN1
Species	Human
Gene ID (Human)	274
Gene ID	274
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-593 of human BIN1 (NP_647601.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00499
Protein Weight	65kDa
Shipping	Ice bag