


## Anti-Bin1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00499

Item Number: CM0006659



### Introduction

Anti-BIN1 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. BIN1 (Amphiphysin II) regulates membrane curvature, endocytosis, and T-tubule formation. Ideal for neuroscience and muscle research.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | BIN1 Rabbit mAb                                                                                                                       |
| Remarks/Alias    | CNM2; AMPH2; AMPHL; SH3P9; BIN1                                                                                                       |
| Species          | Human                                                                                                                                 |
| GeneID (Human)   | 274                                                                                                                                   |
| GeneID           | 274                                                                                                                                   |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 511-593 of human BIN1 (NP_647593.1) |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                              |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | O00499                                                                                                                                |
| Protein Weight   | 65kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |