


Anti-Bin1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00499

Item Number: CM0020810



Introduction

Rabbit polyclonal BIN1 antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human and rat BIN1 (~65 kDa). Ideal for membrane remodeling, endocytosis, and neuromuscular disease research.

[Learn More](#)

Parameter	Value
Product Name	BIN1 Rabbit pAb
Remarks/Alias	CNM2; AMPH2; AMPHL; SH3P9; BIN1
Species	Human
GeneID (Human)	274
GeneID	274
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 230-409 of human BIN1 (NP_647601.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O00499
Protein Weight	65kDa
Shipping	Ice bag