

Anti-Bin3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nqy0

Item Number: CM0018445

Introduction

Anti-BIN3 rabbit polyclonal antibody for Western blot (WB) and ELISA applications. Cross-reacts with human, mouse, and rat BIN3 (~30 kDa). Ideal for cell division and actin cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	BIN3 Rabbit pAb
Remarks/Alias	BIN3
Species	Human
Gene ID (Human)	55909
Gene ID	55909
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-253 of human BIN3 (NP_061158.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQY0
Protein Weight	30kDa
Shipping	Ice bag