

Anti-Biotin Rabbit Monoclonal Antibody For Wb And Elisa

Item Number: CM0000021

Introduction

High-specificity biotin rabbit monoclonal antibody for Western blot (WB) and ELISA. Validated with synthetic peptide immunogen. Wide range predicted cross-reactivity. Ideal for biotin detection in research and IVD workflows.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit mAb
Species	Human
GeneID	CAS:58-85-5
Immunogen	Synthetic peptide Biotin
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Other (Wide Range Predicted)
Shipping	Ice bag