

Anti-Biotinidase (Btd) Polyclonal Antibody For Wb, Ihc-P, Elisa - P43251

Item Number: CM0029710

Introduction

BTD Rabbit pAb targets human biotinidase, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for biotin metabolism research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	BTD Rabbit pAb
Remarks/Alias	BTD
Species	Human
GeneID (Human)	686
GeneID	686
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 284-543 of human BTD (NP_000051.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43251
Protein Weight	61kDa
Shipping	Ice bag