


Anti-Bip/Grp78 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P11021

Item Number: CM0012156



Introduction

Mouse monoclonal antibody targeting human BiP/GRP78, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for endoplasmic reticulum stress, protein folding, and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	BiP/GRP78 Mouse mAb
Remarks/Alias	BIP; GRP78; HEL-S-89n
Species	Human
GeneID (Human)	3309
GeneID	3309
Immunogen	Synthetic peptide A synthesized peptide derived from human BiP/GRP78.
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11021
Protein Weight	72kDa
Shipping	Ice bag