


Anti-Bip/Grp78 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P11021

Item Number: CM0005586



Introduction

High-quality rabbit monoclonal antibody targeting human BiP/GRP78, validated for WB, IF/ICC, IHC-P, and ELISA. Detects endogenous BiP in human, mouse, rat, and pig. Essential for ER stress, chaperone, and viral entry studies.

[Learn More](#)

Parameter	Value
Product Name	BiP/GRP78 Rabbit mAb
Remarks/Alias	BIP; GRP78; HEL-S-89n; BiP/GRP78
Species	Human
Gene ID	3309
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 531-630 of human BiP/GRP78 (NP_005338.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Pig
SWISS	P11021
Protein Weight	72 kDa
Shipping	Ice bag