


Anti-Bip/Grp78 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P11021

Item Number: CM0016258



Introduction

High-quality rabbit polyclonal antibody targeting human BiP/GRP78 (HSPA5). Suitable for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for ER stress and chaperone studies.

[Learn More](#)

Parameter	Value
Product Name	BiP/GRP78 Rabbit pAb
Remarks/Alias	BIP; GRP78; HEL-S-89n; BiP/GRP78
Species	Human
Gene ID (Human)	3309
Gene ID	3309
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 550-650 of human BiP/GRP78 (NP_005338.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11021
Protein Weight	72 kDa
Shipping	Ice bag