


Anti-Bip/Grp78 Rabbit Polyclonal Antibody (Wb, If/Icc, Ihc-P, Elisa) - P11021

Item Number: CM0013569



Introduction

Validated BiP/GRP78 rabbit polyclonal antibody for Western blot (WB), IF/ICC, IHC-P, and ELISA. Recognizes human, mouse, rat BiP, an ER chaperone involved in protein folding and UPR. Ideal for ER stress and cancer research.

[Learn More](#)

Parameter	Value
Product Name	BiP/GRP78 Rabbit pAb
Remarks/Alias	BIP; GRP78; HEL-S-89n; BiP/GRP78
Species	Human
GeneID (Human)	3309
GeneID	3309
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 550-654 of human BiP/GRP78 (NP_005338.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11021
Protein Weight	72 kDa
Shipping	Ice bag