


Anti-Birc7 (Kiap/Livin/ML-Iap) Rabbit Polyclonal Antibody - Q96Ca5

Item Number: CM0012990



Introduction

Rabbit polyclonal antibody against human BIRC7 (Livin, KIAP). Suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for apoptosis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	BIRC7 Rabbit pAb
Remarks/Alias	KIAP; LIVIN; MLIAP; RNF50; ML-IAP; BIRC7
Species	Human
Gene ID (Human)	79444
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human BIRC7 (NP_647478.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96CA5
Protein Weight	33kDa
Shipping	Ice bag